

DETERMINE YOUR APPROPRIATE LENS STRENGTH

THINOPTICS lenses are available in four strength ranges.

Low Range STRENGTH +1.00 ALSO FITS +0.75 to +1.25	Low-Mid Range STRENGTH +1.50 ALSO FITS +1.25 to +1.75	Mid-High Range STRENGTH +2.00 ALSO FITS +1.75 to +2.25	High Range STRENGTH +2.50 ALSO FITS +2.25 to +2.75
--	--	---	---

STEP 1 *Print* this eye test chart. Make sure "Fit to Page" is off. The page should print at actual size. Do not test vision on screen. To make sure this test is accurate, the chart below should be 5 inches wide.

STEP 2 *With your glasses off*, read the lines of copy in the chart below from a distance of about 14 inches.

STEP 3 The first line you have difficulty reading corresponds to the lens strength that is best suited for you.

Printed diagram should be 5" wide.

Ultra-flexible bridge.	+2.75	HIGH RANGE
Instant vision on the go.	+2.50	
Reading glasses in an instant.	+2.25	
NEVER forget your reading glasses again.	+2.00	MID-HIGH RANGE
Full-sized lenses available in four strengths.	+1.75	
Ultra-thin phone case protects your smartphone.	+1.50	
THINOPTICS are insanely flexible and thin reading glasses.	+1.25	LOW-MID RANGE
THINOPTICS featherlight frames come in 6 colors and custom design.	+1.00	
		LOW RANGE



THINOPTICS
ALWAYS WITH YOU

THINOPTICS.COM

Facebook.com/thinoptics Twitter.com/thinoptics

Designed in Palo Alto California from components made in US, Taiwan and China. US and international patents issued and pending. See www.thinoptics.com/patents. Copyright, 2013-2017. THINOPTICS and Always With You are trademarks of THINOPTICS, Inc. LENS 161214